

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
26 August 2004 (26.08.2004)


PCT

(10) International Publication Number
WO 2004/071577 A1

- (51) International Patent Classification⁷: A61N 2/00, 2/12, A61H 23/02
- (21) International Application Number: PCT/US2003/004546
- (22) International Filing Date: 12 February 2003 (12.02.2003)
- (25) Filing Language: English
- (26) Publication Language: English
- (71) Applicant (for all designated States except US): CYMATICS INTERNATIONAL, INC. [US/US]; P.O. BOX 572, PORT JEFFERSON, NY 11777 (US).
- (72) Inventor; and
- (75) Inventor/Applicant (for US only): ARDIZZONE, Vincent [US/US]; CYMATICS INTERNATIONAL, INC., 9 STUARTS CT, PORT JEFFERSON, NY 11777 (US).
- (81) Designated States (national): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.
- (84) Designated States (regional): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

[Continued on next page]

(54) Title: MAGNETO-CYMATIC THERAPEUTIC OSCILLATING WAND



(57) Abstract: A magneto-cymatic therapeutic device providing simultaneously magnetic therapy with cymatic therapy to an anatomical area to be treated, delivering the cymatic frequency via mechanical and/or magnetic interaction. A transducer coupled to a magnet converts audio frequency electrical signals into mechanical vibrations, within an applicator-wand of a portable magneto-cymatic therapeutic device. To treat an anatomical area the applicator-wand may be rotated, or moved through a particular pattern. The application surface of the applicator-wand may be contoured to allow maximum contact with the anatomical body part to be treated. The magnetic field extending out from the application surface also interacts with tissue in the body and also delivers cymatic vibrations without physical contact between the device and the body.



Published:

— with international search report

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.